

RI-MAGIC LED LIGHT - trolley - 6250



Product Code	49051
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396168275

Description

RI-MAGIC® HP LED - trolley

High performance medical LED light for examination and procedures.

- high performance LED, 6,000 °K
- lifetime of the lamp up to 50,000 hours
- significantly reduced heating and energy loss
- focusing device at the front of the light head
- highly flexible, long special arm can be bent in all directions
- rheostat handle for infinitely variable regulation of light intensity.

Manual in: GB, FR, IT, ES, DE, RU.

Sturdy mobile base made of fiber glass and chrome plated pole
Low friction antistatic castors with brake
Ø 53 cm mobile base

Made in Germany.

Technical Specifications

Working distance: 200 mm

Lux max : 55,000 lux

Circle of light max Ø: 90 mm

Circle of light min Ø: 20 mm

Working distance: 400 mm

Lux max : 45,000 lux

Circle of light max Ø: 155 mm

Circle of light min Ø: 25 mm

Working distance: 600 mm

Lux max : 30,000 lux

Circle of light max Ø: 215 mm

Circle of light min Ø: 30 mm

Working distance: 800 mm

Lux max : 20,000 lux

Circle of light max Ø: 270 mm

Circle of light min Ø: 35 mm

Working distance: 1,000 mm

Lux max : 15,000 lux

Circle of light max Ø: 340 mm

Circle of light min Ø: 40 mm

Operating Voltage: 100-240 V • 50-60 Hz