



Product Code	49053
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396195813

Description

RI-MAGIC® HP LED - trolley - 6252 with power pack fixed on stand

High performance medical LED light for examination and procedures.

- high performance LED, 6,000 °K

- lifetime of the lamp up to 50,000 hours

- significantly reduced heating and energy loss

- focusing device at the front of the light head

- highly flexible, long special arm can be bent in all directions

- rheostat handle for infinitely variable regulation of light intensity.

Sturdy mobile base made of fiber glass and chrome plated pole.

Manual in: GB, FR, IT, ES, DE, RU.

Made in Germany.

Low friction antistatic castors with brake

Ø 53 cm mobile base

Technical Specifications

Working distance: 200 mm

Lux max : 55,000 lux

Circle of light max Ø: 90 mm

Circle of light min Ø: 20 mm

Working distance: 400 mm

Lux max : 45,000 lux

Circle of light max Ø: 155 mm

Circle of light min Ø: 25 mm

Working distance: 600 mm

Lux max : 30,000 lux

Circle of light max Ø: 215 mm

Circle of light min Ø: 30 mm

Working distance: 800 mm

Lux max : 20,000 lux

Circle of light max Ø: 270 mm

Circle of light min Ø: 35 mm

Working distance: 1,000 mm

Lux max : 15,000 lux

Circle of light max Ø: 340 mm

Circle of light min Ø: 40 mm

Operating Voltage: 100-240 V • 50-60 Hz