



Product Code	34908
Unit of sale	50 bags of 6 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2183488
CND	H900102
PARAF	977814957
GTIN	30707387804641

## Description

### 3M STERI-STRIP™ - SKIN CLOSURE STRIPS

#### REINFORCED STERI-STRIP™

With polyester reinforcing filaments for higher tensile strength, reduced infection and scarring risk and better cosmetic outcomes.

- STERI-STRIP™ - 6 x 38 mm
- Reinforced
- Sterile
- 3M Code: R1542
- Minimum order: 50 bags of 6 strips

Adhesive skin closures made of non-woven backing coated with a pressure-sensitive hypoallergenic adhesive. Sterile, breathable and comfortable to wear.