



Product Code	49133
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2723948
CND	Z12010701
EAN/UPC	8436562860431

Description

MS LED FLEX PLUS LIGHT - trolley

Minor surgery luminaire with LED light source. Flexible or articulated arm with great precision allow regulation of the light in any position. High quality of materials makes it especially suitable for intensive use in surgery and examinations.

Exceptional quality of LED light

Powerful LED light source offers exceptional light quality and a lifespan of more than 50,000 hours. The high performance of the colour Ra 97% depicts the exceptional quality of white light for an excellent distinction of fabrics and reproduction of all visible colours, naturally and precisely.

Ergonomic design

The circular shape of the head and its handle in the lower part facilitate its exact positioning, clean and soft shapes which facilitate hygiene and cleanliness.

Dimmable, adjustable light intensity

Made in Spain.

Technical Specifications

Light source LED 12W
 Colour rendering > 97%
 Lux at 50 centimeters 65.000 lux
 Lamp lifetime 50.000 hours
 Adjustable light intensity
 Colour temperature 4.500 °K
 Arm length 70 cm
 Colour White