

TEGADERM 3M I.V. FIXING STRIPS 7x8.5 cm



Product Code	34945
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2265709
CND	M04010202
EAN/UPC	8711428084665
PARAF	977815164

Description

3M TEGADERM™ IV - STERILE - 7 x 8.5 cm - for peripheral lines
3M code: 1633

TEGADERM™ IV
A cost-effective solution for peripheral and central vascular lines with strong hold and gentle removal.

A complete range of sterile, waterproof, highly-breathable transparent film dressings to protect IV sites from water, virus and bacteria. With antifriction backing and reinforced non-woven borders for longer wearing time and minimized catheter dislodgement risk.