



|               |                   |
|---------------|-------------------|
| Product Code  | 21495             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

EAN/UPC 8023279214956

## Description

TROUSERS - UNISEX - 60% COTTON - 40% POLYESTER

- Size: XL
- Colour: white
- UNISEX
- Elasticated waist
- 2 front and 1 back pockets

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.  
Stock in 5 colours, other 22 colours available in 60 days.

## Technical Specifications

Composition: 60% cotton - 40% polyester  
Weight: 200 g/m2 (white), 185 g/m2 (coloured)  
Yarn count: warp: Ne 20 ± 2 weft: Ne 18 ± 2  
Construction: warp: 40 ± 3 weft: 19 ± 3  
Breaking strength - warp: 100 kg ±3%  
Breaking strength - weft: 39 kg ±3% (white) 49 kg ±3% (coloured)  
Fastness to washing: 4  
Treatments and finishing: crease resist  
Dimensional stability - washing at 90°C: 3%