

TRANSPORE 3M WHITE TAPE 13 mm x 9.14 m



Product Code	34950
Unit of sale	box of 24 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2192657
CND	M050101
PARAF	977815188
GTIN	04064035092665

Description

3M™ Transpore™ White Surgical Tape

- 3M code: 1534-0
- Dimensions: 13 mm x 9.14 m

Gentle tape with strong adhesion and easy, bi-directional tear for dressing applications.

Hypoallergenic, breathable, easy bi-directional tear. Good initial and long-term adhesion (72 hours), easy to handle with gloves, water-resistant.

Suggested Applications:

- securing dressings to damp or dry skin
- taping fragile, at-risk skin
- anchoring tubes and devices
- when repeated taping is needed
- when tape width must be customized