

YTON NURSE STETHOSCOPE - Y dark blue



Product Code	49501
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2513120
CND	C9005
EAN/UPC	8023279495010
PARAF	952100269
GTIN	08023279495010
GMDN	13755

Description

YTON STETHOSCOPES - LATEX FREE

- Yton single head - dark blue

Solid aluminium chestpiece Ø 45 mm, bell Ø 30 mm (for double head), with colour coordinated plastic retaining ring.
"Duracote" finished.

Available in four different colours: dark green, dark blue, red and aqua blue.

Manual and box in 11 languages: GB, FR, IT, DE, ES, PT, GR, PL, SE, HU, Arabic.