

TRANSPORE 3M WHITE TAPE 25 mm x 9.14 m



Product Code	34951
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2192658
CND	M050101
EAN/UPC	4064035094829
PARAF	977815190
GTIN	04064035094829

Description

3M™ Transpore™ White Surgical Tape

- 3M code: 1534-1
- Dimensions: 25 mm x 9.14 m

Gentle tape with strong adhesion and easy, bi-directional tear for dressing applications.

Hypoallergenic, breathable, easy bi-directional tear. Good initial and long-term adhesion (72 hours), easy to handle with gloves, water-resistant.

Suggested Applications:

- securing dressings to damp or dry skin
- taping fragile, at-risk skin
- anchoring tubes and devices
- when repeated taping is needed
- when tape width must be customized