

YTON DUAL HEAD STETHOSCOPE - Y dark green



Product Code	49510
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2513127
CND	C9005
EAN/UPC	8023279495102
GTIN	08023279495102
GMDN	13755

Description

YTON STETHOSCOPES - LATEX FREE

- Yton double head - dark green

Solid aluminium chestpiece Ø 45 mm, bell Ø 30 mm (for double head), with colour coordinated plastic retaining ring.
"Duracote" finished.

Available in four different colours: dark green, dark blue, red and aqua blue.

Manual and box in 11 languages: GB, FR, IT, DE, ES, PT, GR, PL, SE, HU, Arabic.