

TEGADERM 3M FOAM 5x5 cm - non adhesive



| | |
|---------------|----------------|
| Product Code | 34953 |
| Unit of sale | box of 10 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|------------|----------------|
| RDM (NSIS) | 2671221 |
| CND | M04040601 |
| GTIN | 00707387807184 |

Description

3M™ TEGADERM™ HP FOAM DRESSINGS

3M code: 90600

Product: Tegaderm™ HP foam

Size: 5x5 cm

Minimum order: Box of 10

Type: Non-adhesive dressing

3M™ TEGADERM™ HIGH PERFORMANCE FOAM

Sterile, multi-layer foam dressing that combines high absorbency with high breathability to reduce the risk of wound maceration.

Well suited for changing levels of exudate to maintain moisture balance for optimal wound healing.

Highly breathable and waterproof, also available in heel/elbow shape with adhesive border and simplified application system.