

TEGADERM 3M FOAM 8.8x8.8 cm - non adhesive



Product Code	34955
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B
RDM (NSIS)	2671226
CND	M040203
EAN/UPC	8023279349559
PARAF	977815226

Description

3M™ TEGADERM™ HP FOAM DRESSINGS

3M code: 90604
Product: Tegaderm™ HP foam
Size: 8.8x8.8 cm
Minimum order: Box of 10
Type: Non-adhesive dressing
Indication: For tracheostomy

3M™ TEGADERM™ HIGH PERFORMANCE FOAM

Sterile, multi-layer foam dressing that combines high absorbency with high breathability to reduce the risk of wound maceration. Well suited for changing levels of exudate to maintain moisture balance for optimal wound healing. Highly breathable and waterproof, also available in heel/elbow shape with adhesive border and simplified application system.