



Product Code	34959
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	III

RDM (NSIS)	2395613
CND	M04010299
GTIN	30707387782178

Description

3M™ TEGADERM™ CHG IV SECUREMENT DRESSING

3M code: 1657R

Product: Tegaderm™ CHG - 8.5x11.5 cm

Minimum order: box of 25

Integrated I.V. dressing combines infection reduction with site visibility, catheter securement, and breathability. Proven to reduce CRBSI and vascular catheter colonization. Immediate and continuous antimicrobial protection through integrated chlorhexidine gluconate gel pad, without requiring moisture to activate it. Designed to minimize catheter movement, it can be worn up to 7 days. Waterproof and sterile barrier to external contaminants (liquids, bacteria, and viruses).