

TEGADERM 3M PICC/CVC + CHG IV ADVANCED

8.5x11.5 cm



Product Code	34963
Unit of sale	box of 20 pcs
Minimum order	1
Type	Kit of medical devices
Class	

RDM (NSIS)	139516
EAN/UPC	0707387769482
GTIN	50707387769487

Description

3M™ TEGADERM™ PICC/CVC IV ADVANCED SECUREMENT DRESSING

- 3M code: 1877R-2100
- Tegaderm™ PICC/CVC + CHG - 8.5x11.5 cm
- box of 20

3M™ TEGADERM™ PICC/CVC advanced IV fixation dressing

Waterproof transparent film which provides a barrier against external contaminants (liquids, bacteria and viruses).

The pre-formed plastic fastener is integrated on a breathable base with a silicone adhesive, that can be repositioned and easily removed from the skin avoiding needle injuries.

They allow for continuous monitoring of the IV insertion site.

Designed to stay in place for up to 7 days.

Suitable for fixation of peripherally inserted central catheter (PICC) and short-term central venous catheters (CVC).