

TEGADERM 3M ALGINATE 10x10 cm



Product Code	34964
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2056606
CND	M04040201
EAN/UPC	4054596862081
PARAF	977815240
GMDN	43186

Description

TEGADERM™ ALGINATE
High integrity, gel forming alginate dressing for moderate to heavily exuding wounds. Compact and comfortable.

- 3M code: 90112
- TEGADERM™ ALGINATE
- 10 x 10 cm