

CUROS DISINFECTING CAPS - needle-free connector



Product Code	34980
Unit of sale	box of 270 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2316217
CND	A0706
EAN/UPC	4054596425446

Description

3M™ CUROS™ DISINFECTING CAPS

- 3M code: CFF1-270R
- CuroS™ caps
- Connector: needle-free
- Sale formats: single
- box of 270

3M™ CUROS™ DISINFECTING CAPS

3M™ CuroS caps disinfect in 1 minute and protect needle connectors and Luer to help reducing the risk of contaminants from entering the catheter post-insertion, for up to 7 days if not removed.
Caps twist on easily and stay securely in place, preventing contained alcohol (70% isopropyl alcohol) from spreading within patient's circulatory system.
For single use only.
Suitable for needle-free connectors (34980), male Luer (34981) and open female Luer (34982).