

HEINE CUFF 1 tube - child - 2



Product Code	49801
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
EAN/UPC	4053755136940
PARAF	987475199
GMDN	34978

Description

HEINE ORIGINAL CUFFS 1 TUBE

- 1 TUBE
- Child - 2
- Heine code: M-000-09-702
- Colour: grey
- Size cm: 10x34
- Circumference cm: 13-20

High quality single tube cuffs with velcro fastener. Perfectly balanced for comfort and ease of use. Compatible with all Heine sphygmomanometers. Water resistant, washable material. Latex free.

Compatibility: All Heine 1 tube