

HEINE CUFF 1 tube - adult small - 3



Product Code	49802
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

EAN/UPC	4053755136964
PARAF	987475187
GMDN	34978

Description

HEINE ORIGINAL CUFFS 1 TUBE

- 1 TUBE
- Adult small - 3
- Heine code: M-000-09-704
- Colour: grey
- Size cm: 14x53
- Circumference cm: 20-29

High quality single tube cuffs with velcro fastener. Perfectly balanced for comfort and ease of use. Compatible with all Heine sphygmomanometers. Water resistant, washable material. Latex free.

Compatibility: All Heine 1 tube