

HEINE CUFF 1 tube - adult - 4



Product Code	49803
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
EAN/UPC	4053755136759
PARAF	987475163
GMDN	34978

Description

HEINE ORIGINAL CUFFS 1 TUBE

- Adult - 4
- Heine code: M-000-09-619
- Colour: grey
- Size cm: 14x58
- Circumference cm: 29-41

High quality single tube cuffs with velcro fastener. Perfectly balanced for comfort and ease of use. Compatible with all Heine sphygmomanometers. Water resistant, washable material. Latex free.

Compatibility: All Heine 1 tube