

CUROS TIPS DISINFECTING CAPS - male Luer connector - 40 strips of 5



Product Code	34981
Unit of sale	box of 200 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2316278
CND	A0706
GTIN	30707387773800

Description

3M™ CUROS™ DISINFECTING CAPS

- 3M code: CM5-200R
- Curostips™ caps
- Connector: male Luer
- Sale formats: 40 strips of 5
- box of 200

3M™ CUROS™ DISINFECTING CAPS

3M™ Curostips caps disinfect in 1 minute and protect needle connectors and Luer to help reducing the risk of contaminants from entering the catheter post-insertion, for up to 7 days if not removed.

Caps twist on easily and stay securely in place, preventing contained alcohol (70% isopropyl alcohol) from spreading within patient's circulatory system.

For single use only.

Suitable for needle-free connectors (34980), male Luer (34981) and open female Luer (34982).