

RIESTER CUFF 1 TUBE - thigh



Product Code	49814
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2808683
CND	Z12040185
EAN/UPC	4045396001022
PARAF	987475720

Description

RIESTER ORIGINAL CUFF - LATEX FREE

- 1 TUBE
- Thigh velcro cuff
- Riester code: 125
- Colour: black
- Size cm: 70x22
- Circumference cm: 42 - 50

Compatibility:
All Riester single tube manual aneroid sphygmomanometers