

RIESTER ONE PIECE CUFF 1 TUBE - child



Product Code	49816
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	1567166
CND	C9006
EAN/UPC	4045396168183

Description

ONE PIECE RIESTER ORIGINAL CUFF - LATEX FREE

- 1 TUBE - length 50 cm
- Child one-piece cuff
- Riester code: 154
- Colour: dark green
- Size cm: 34x10
- Circumference cm: 13 - 17

Compatibility:
All Riester single tube manual aneroid sphygmomanometers

They prevent cross contaminations. Comfortable soft surface, latex free.
Washable at 60°C and disinfectable in liquid, reusable up to 10,000 times. ISO 105 E01 tested.