


RIESTER CUFF 2 TUBES - adult small



Product Code	49820
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396000988

Description

RIESTER ORIGINAL CUFF - LATEX FREE

- 2 TUBES
- Adult small velcro cuff
- Riester code: 121
- Colour: black
- Size cm: 42x13
- Circumference cm: 17 - 24

Compatibility:

All Riester twin tube sphygmomanometers including BigBen