



| | |
|---------------|----------------|
| Product Code | 49821 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|---------|---------------|
| EAN/UPC | 4045396001077 |
|---------|---------------|

Description

RIESTER ORIGINAL CUFF - LATEX FREE

- 2 TUBES
- Child velcro cuff
- Riester code: 130
- Colour: light blue
- Size cm: 35.5x10
- Circumference cm: 13 - 17

Compatibility:

All Riester twin tube sphygmomanometers including BigBen