

RIESTER CUFF 2 TUBES - adult large



Product Code	49823
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

CND	V0399
EAN/UPC	4045396001008
PARAF	987475769

Description

RIESTER ORIGINAL CUFF - LATEX FREE

- 2 TUBES
- Adult large velcro cuff
- Riester code: 123
- Colour: black
- Size cm: 70x15
- Circumference cm: 34 - 44

Compatibility:
All Riester twin tube sphygmomanometers including BigBen