

# RIESTER CUFF 2 TUBES - thigh



Product Code	49824
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

CND	V0399
EAN/UPC	4045396001039
PARAF	987475732

## Description

### RIESTER ORIGINAL CUFF - LATEX FREE

- 2 TUBES
- Thigh velcro cuff
- Riester code: 126
- Colour: black
- Size cm: 70x22
- Circumference cm: 42 - 50

#### Compatibility:

All Riester twin tube sphygmomanometers including BigBen