

RIESTER ONE PIECE CUFF 2 TUBES - adult



Product Code	49827
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396168152
PARAF	987475744

Description

ONE PIECE RIESTER ORIGINAL CUFF - LATEX FREE

- 2 TUBES - length 50 cm
- Adult one-piece cuff
- Riester code: 151
- Colour: dark green
- Size cm: 54.5x14.5
- Circumference cm: 24 - 32

Compatibility:

All Riester twin tube sphygmomanometers including BigBen

They prevent cross contaminations. Comfortable soft surface, latex free.
Washable at 60°C and disinfectable in liquid, reusable up to 10,000 times. ISO 105 E01 tested.