

# ADAPTIC TOUCH 3M NON-ADHERING SILICONE DRESSING 7.6x5 cm



Product Code	34985
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2703326
CND	M04040702
GTIN	04054596951594

## Description

### 3M™ ADAPTIC™ DRESSINGS

3M™ code: TCH501  
Product: Adaptic™ Touch - silicone  
Size: 7.6x5 cm  
Minimum order: Box of 10

3M™ Adaptic™ are primary wound contact dressings. Pore size is designed to minimize the pooling of fluid at the wound site and enables fluid transfer to a secondary dressing, reducing the risk of maceration. They minimise adherence and trauma associated with dressing changes, protecting newly-formed tissue and reducing pain. Ideal for low to high levels of exudate (ulcers, first- and second-degree burns, abrasions...).

3M™ Adaptic™ Touch Dressing is made of soft tack silicone, also suitable for use, under medical supervision, with negative pressure wound therapy (NPWT) for the protection of fragile wound structures.