

BP2A BLOOD PRESSURE MONITOR



Product Code	49865
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
RDM (NSIS)	2733178
CND	Z1203020501
EAN/UPC	6934440102304
PARAF	950975918
GMDN	45617

Description

BP2A Bluetooth Blood Pressure Monitor

Intelligent automatic upper arm blood pressure monitor designed for daily use. The all-in-one design makes it easy to carry. Clinically validated: combination of Intelligent pressurization technology and the brand-new Lycra high elastic cuff design ensure the most accurate measurement results.

Easy to read and use: one-click measurement and OLED screen shows real-time monitoring.

Reliability: detects body movement and cuff wearing to avoid invalid results.

Bluetooth sync allows monitoring, sharing, and saving results on ViHealth App, which supports unlimited storage and multiuser management.

Supplied with integrated cuff (22-42 cm), built-in rechargeable battery, charging cable, quick guide, and multilingual manual in: GB, FR, IT, ES, DE.

Technical Specifications

Measuring method: oscillometric method
Pressure measur. range: 0-300 mmHg; ± 3 mmHg or 2% whichever is greater
Pulse rate range: 40 to 200 bpm; ± 2 bpm
Storage: 50 records
Charge input: USB Type-C, DC 5V
Power supply: Li-ion rechargeable battery
Battery life: 500 measurements
Charging time: 2 hours
Bluetooth: Bluetooth 4.0 BLE
Size - Weight: 135x45xh 20 mm - 240 g