

X-CHECK PRO AUTOMATIC BLOOD PRESSURE MONITOR with Bluetooth



| | |
|---------------|----------------|
| Product Code | 49872 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2761614 |
| CND | Z1203020501 |
| EAN/UPC | 6972040116185 |

Description

X-CHECK PRO AUTOMATIC B.P. MONITOR with Bluetooth

Automatic blood pressure monitor which combines smart technology with medical-grade accuracy to make home monitoring easier and more reliable. It features a large 4.2" LCD display with a clear 4-colour backlight system. The upgraded algorithm and intelligent inflation technology ensure precise readings, while the irregular heartbeat (IHB) detection helps identify abnormal heart rhythms during each measurement. A cuff wearing detection system confirms when the cuff is properly positioned with an "OK" indicator, and movement detection ensures that measurements are taken under the right conditions.

- one-button operation

- 2x120 memories with date and time

- average reading of last 3 measurements

It pairs with the free Tuya Smart App, to track and manage health data.

Supplied with adult soft nylon cuff 22-42 cm, pouch and multilingual manual in GB, FR, IT, ES, PT, DE, NL, SE, PL, HU, RO, GR.

Technical Specifications

Measuring method: oscillometric

Display size: 80x68.9 mm

Measuring range pressure: SYS 57-255 mmHg / DIA 25-195 mmHg

Measuring range pulse: 40-199 bpm

Memory: 2x120 memories

Power supply: 4xAAA batteries or type-C USB cable (both not included)

Power off: automatic switch off after 60 seconds

Display: digital LCD display, 4-colour backlight

Size - weight: 108x139x62 mm - 316 g (without batteries)