

OMRON EVOLV ALL IN ONE ARM BLOOD PRESSURE MONITOR



Product Code	49900
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1496402
CND	Z1203020501
EAN/UPC	4015672110465
PARAF	971992413

Description

The new, easy to use, all in one upper arm blood pressure monitor.

Take accurate readings in any position around the upper arm and track progress via smartphone.

- All in one upper arm blood pressure monitor: no tubes, no wires, simply accurate
- Intellisense Technology: inflates the cuff to the ideal level
- Intelli Wrap cuff: for an easy and accurate reading
- Clinically Validated: proven reliability and accuracy under vigorous testing procedures by major health organisations.
- Cuff Wrap guide: indicates if the cuff is wrapped correctly, not too loose or too tight
- OMRON connect App: get connected, access anytime, see progress. Mobile application for iPhone / Android™.

Supplied with battery set, storage case, manuals.
Compatible with Alexa.

Multilingual box and manual: GB, FR, IT, ES, DE, NL, RU, TR, Arabic.

With standard cuffs, 1-in-3 people measure their blood pressure incorrectly. EVOLV's Intelli Wrap Cuff Technology ensures 360° accuracy in any position around your upper arm.

EVOLV can sync your readings to your smartphone to track results, monitor progress, and get a deeper understanding of your health.
See smartphone compatibility list on omronconnect.com/emea/en_gb/

Technical Specifications

Cuff Type: Intelli Wrap Cuff 22-42 cm
Memory: 100 Readings
Cuff Wrap Guide: Yes
Irregular Heartbeat Detection: Yes
Body Movement Detection: Yes
Validation: Clinical Validation, Pregnancy Validation
Intellisense Technology: Yes
Stroke Prevention: Yes
Size: 12x17.3x12.5 cm