



Product Code	49903
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2156655
CND	Z1203020501
EAN/UPC	4015672113046

Description

OMRON M2 INTELLI IT HEM-7143T-EBK

With cuff 22-32 cm. Unique Intellisense™ Technology and Irregular Heartbeat Detection. One-touch design and large display.

Clinically validated (To verify validations, please see www.dablededucational.org)

- Intellisense technology
- Irregular heartbeats detector
- Cuff wrap guide
- Cuff: 22-32 cm
- Memory measurem. x user: 30
- Body movement detection
- Hypertension
- App/Bluetooth Connection

Technical Specifications

- Measuring method: oscillometric
- Measuring range - Pressure: 0 mmHg to 299 mmHg \pm 3 mmHg
- Measuring range - Pulse: 40 to 180 min \pm 5% of display reading
- Memory (with date and time): 30
- Inflation: fuzzy logic controlled by electric pump
- Deflation: automatic pressure release valve
- Rating: DC 6V 4W
- Power source: 4 "AA" batteries 1.5 V
- Battery life (N° of measurements): 1000
- Protection against electric shock: internally powered ME equipment (when using batteries)
- Weight (console+cuff) without batteries: 250+130 g
- Size: 103x140xh 82 mm
- Cuff materials: nylon and polyester
- Measurable circumference (cm): 22-32 (arm)