



non inclus
not included
non inclus



Product Code	49904
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2817605
CND	Z1203020501
EAN/UPC	4015672114289

Description

OMRON M2+ CONNECT HEM-7188T1-LE

With cuff 22-42 cm. Unique Intellisense™ Technology and Irregular Heartbeat Detection. One-touch design and large display.

Clinically validated (To verify validations, please see www.dablededucational.org) Intellisense technology Irregular heartbeats detector Cuff wrap guide Cuff model: Large Cuff - Size: 22-42 cm Memory measurement. x user: 30 Body movement detection Easy high blood pressure colour indicator (Signal in case the values exceed the E.H.S. standards for domiciliary auto-measurements (85 systolic pressure and/or 135 diastolic pressure)) Hypertension App/Bluetooth Connection

The Intellisense™ technology ensures that the right inflation pressure is automatically detected. It detects irregular pulses and automatically decides if the measurement result is reliable or needs to be repeated. Only accurate results are displayed.

Technical Specifications

Measuring method: oscillometric

Measuring range: Pressure: 0 mmHg to 299 mmHg ± 3 mmHg Pulse: 40 to 180 min ± 5% of display reading

Memory (with date and time): 30 Inflation: fuzzy logic controlled by electric pump Deflation: automatic pressure release valve Rating: DC 6V 4W Power source: 4 "AA" batteries 1.5 V Battery life (N° of measurements): 900 Protection against electric shock: internally powered ME equipment (when using batteries) Weight (console+cuff) without batteries: 250+130 g Size: 103x140xh 82 mm Cuff materials: nylon and polyester Measurable circumference (cm): 22-42 (arm)

Standard accessories

Standard languages: GB, FR, IT, ES, DE, NL, RU, TR, Arabic Other languages: available on Omron website