



Product Code	49905
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2556510
CND	Z1203020501
EAN/UPC	4015672111837
PARAF	981516103

Description

OMRON M3 HEM-7154-E with Easy cuff 22-42 cm.

- Clinically validated: to verify validations, please see www.dablededucational.org
- Intellisense Technology
- Irregular heartbeat detector
- Cuff wrap guide
- Intelli Wrap Technology: NO
- Cuff: Easy cuff 22-42 cm
- Memory measurement: 60 x 2 users + guest mode
- Body movement detection
- Easy high blood pressure colour indicator: signal in case the values exceed the E.H.S. standards for domiciliary auto-measurements (85 systolic pressure and/or 135 diastolic pressure)
- Average last 3 measurements (Average of 3 last measurements made within 10 min)

The Intellisense™ technology ensures that the right inflation pressure is automatically detected.

It detects irregular pulses and automatically decides if the measurement result is reliable or needs to be repeated.

Only accurate results are displayed.

Multilingual box and manual: GB, FR, IT, ES, DE, NL, RU, TR, Arabic.

Technical Specifications

Measuring method: oscillometric

Measuring range:

- Pressure: 0 mmHg to 299 mmHg ± 3 mmHg

- Pulse: 40 to 180 min ± 5% of display reading

Memory (with date and time): 60x2 users

Inflation: fuzzy logic controlled by electric pump

Deflation: automatic pressure release valve

Rating: 6VCC 4W

Power source: 4 "AA" batteries 1.5 V

Battery life (N° of measurements): 1000

Protection against electric shock: internally powered ME equipment (when using batteries)

Weight (console+cuff) without batteries: 310+170 g

Size: 105x153xh 87 mm

Cuff materials: nylon and polyester

Measurable circumference (cm): 22-42 (arm)