



Product Code	49908
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2556625
CND	Z1203020501
EAN/UPC	4015672111776

Description

OMRON M4 INTELLI IT HEM-7155T-EBK

With Intelli Wrap Cuff 22-42 cm.
OMRON App to monitor hypertension at home.

The Intellisense™ technology ensures that the right inflation pressure is automatically detected. It detects irregular pulses and automatically decides if the measurement result is reliable or needs to be repeated. Only accurate results are displayed.

- Clinically validated (to verify validations, please see www.dableducational.org)
- Intellisense technology
- Irregular heartbeats detector
- Cuff wrap guide
- Intelli Wrap Technology (360° accuracy: accurate measurements from any position around the upper arm)
- Pre-formed cuff: Intelli wrap 22-42 cm
- Memory measurem. x user: 60x2
- Guest mode
- Body movement detection
- Easy high blood pressure colour indicator. Signal in case the values exceed the E.H.S. standards for domiciliary auto-measurements (85 systolic pressure and/or 135 diastolic pressure)
- Hypertension
- Average last 3 measurements (Average of 3 last measurements made within 10 min)
- Measurements comparison versus prior readings
- App/Bluetooth Connection

Technical Specifications

- Measuring method: oscillometric
- Measuring range:
Pressure: 0 mmHg to 299 mmHg ± 3 mmHg
Pulse: 40 to 180 min ± 5% of display reading
- Memory (with date and time): 60x2 users
- Inflation: fuzzy logic controlled by electric pump
- Deflation: automatic pressure release valve
- Rating: 6VCC 4W
- Power source: 4 "AA " batteries 1.5 V
- Battery life (N° of measurements): 1000
- Protection against electric shock: internally powered ME equipment (when using batteries)
- Weight (console+cuff) without batteries: 340+163 g
- Size: 105x152xh 85 mm
- Cuff materials: nylon and polyester
- Measurable circumference (cm): 22-42 (arm)