



Product Code	49911
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2576435
CND	Z1203020501
EAN/UPC	4015672113978
PARAF	988050086

Description

OMRON M6 COMFORT - AFIB - HEM-7380-E

With Intelli Wrap Cuff 22-42 cm.
It detects the possibility of Atrial fibrillation, enabling home monitoring. It also takes 3 consecutive readings at 30 second intervals, then displays the average to give more accurate indication of blood pressure.

The Intellisense™ technology ensures that the right inflation pressure is automatically detected. It detects irregular pulses and automatically decides if the measurement result is reliable or needs to be repeated. Only accurate results are displayed.

- Clinically validated (To verify validations, please see www.dablededucational.org)
- Intellisense technology
- Irregular heartbeats detector (AFib)
- Cuff wrap guide
- Intelli Wrap Technology (360° accuracy: accurate measurements from any position around the upper arm)
- Pre-formed cuff: model Intelli wrap / size (cm) 22-42 (Pre-formed cuff, excluding M2+, M2 Intelli IT, M3)
- Memory measurem. x user 100x2
- Guest mode
- Body movement detection
- Easy high blood pressure colour indicator (Signal in case the values exceed the E.H.S. standards for domiciliary auto-measurements (85 systolic pressure and/or 135 diastolic pressure))
- Hypertension
- Average last 3 measurements (Average of 3 last measurements made within 10 min)
- AFib detection in each measurement
- Measurements comparison versus prior readings

Technical Specifications

- Measuring method: oscillometric
- Measuring range:
 - Pressure: 0 mmHg to 299 mmHg ± 3 mmHg
 - Pulse: 40 to 180 min ± 5% of display reading
- Memory (with date and time): 100x2 users
- Inflation: fuzzy logic controlled by electric pump
- Deflation: automatic pressure release valve
- Rating: 6VCC 4W
- Power source: 4 "AA " batteries 1.5 V
- Battery life (N° of measurements): 1000
- Protection against electric shock: internally powered ME equipment (when using batteries)
- Weight (console+cuff) without batteries: 460+163 g
- Size: 191x120xh 85 mm
- Cuff materials: nylon and polyester
- Measurable circumference (cm): 22-42 (arm)