



Product Code	49915
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2556619
CND	Z1203020501
EAN/UPC	4015672111721
PARAF	981516178

Description

OMRON M7 INTELLI IT - HEM-7361T-EBK with Intelli Wrap Cuff 22-42 cm.

- Clinically validated: to verify validations, please see www.dablededucational.org
- Intellisense Technology
- Irregular heartbeat detector: AFib
- Cuff wrap guide
- Intelli Wrap Technology (360° accuracy: accurate measurements from any position around the upper arm)
- Cuff: Intelli wrap cuff 22-42 cm
- Memory measurement: 100 x 2 users + guest mode
- Body movement detection
- Easy high blood pressure colour indicator: signal in case the values exceed the E.H.S. standards for domiciliary auto-measurements (85 systolic pressure and/or 135 diastolic pressure)
- Hypertension
- Average last 3 measurements (Average of 3 last measurements made within 10 min)
- 3 consecutive measurements in AFib mode
- Measurements comparison versus prior readings
- App/Bluetooth Connection

The Intellisense™ technology ensures that the right inflation pressure is automatically detected.

It detects irregular pulses and automatically decides if the measurement result is reliable or needs to be repeated.

Only accurate results are displayed.

The connect OMRON App allows to monitor parameters, create personalized programs, download and share data. AFib (Atrial Fibrillation) detection function.

Multilingual box and manual: GB, FR, IT, ES, DE, NL, RU, TR, Arabic.

Technical Specifications

Measuring method: oscillometric

Measuring range:

- Pressure: 0 mmHg to 299 mmHg ± 3 mmHg

- Pulse: 40 to 180 min ± 5% of display reading

Memory (with date and time): 100x2 users

Inflation: fuzzy logic controlled by electric pump

Deflation: automatic pressure release valve

Rating: 6VCC 4W

Power source: 4 "AA" batteries 1.5 V

Battery life (N° of measurements): 1000

Protection against electric shock: internally powered ME equipment (when using batteries)

Weight (console+cuff) without batteries: 460+163 g

Size: 191x120xh 85 mm

Cuff materials: nylon and polyester

Measurable circumference (cm): 22-42 (arm)