

OMRON M7 INTELLI IT AFIB DIGITAL B.P.M. HEM-7380T1-EBK



Product Code	49916
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2576561
CND	Z1203020501
EAN/UPC	4015672113923
PARAF	988050098

Description

OMRON M7 INTELLI IT - AFIB HEM-7380T1 -EBK
With Intelli Wrap Cuff 22-42 cm.

The connection to OMRON App allows to monitor parameters, create personalized programs, download and share data. AFib (Atrial Fibrillation) detection function.

The Intellisense™ technology ensures that the right inflation pressure is automatically detected. It detects irregular pulses and automatically decides if the measurement result is reliable or needs to be repeated. Only accurate results are displayed.

Clinically validated (To verify validations, please see www.dablededucational.org)

Intellisense technology

Irregular heartbeats detector (AFib)

Cuff wrap guide

Intelli Wrap Technology (360° accuracy: accurate measurements from any position around the upper arm)

Pre-formed cuff: model Intelli wrap / size (cm) 22-42 (Pre-formed cuff, excluding M2+, M2 Intelli IT, M3)

Memory measur. x user 100x2

Guest mode

Body movement detection

Easy high blood pressure colour indicator (Signal in case the values exceed the E.H.S. standards for domiciliary auto-measurements (85 systolic pressure and/or 135 diastolic pressure))

Hypertension

Average last 3 measurements (Average of 3 last measurements made within 10 min)

AFib detection in each measurement

Measurements comparison versus prior readings

App/Bluetooth Connection

Technical Specifications

Measuring method: oscillometric

Measuring range:

- Pressure: 0 mmHg to 299 mmHg \pm 3 mmHg

- Pulse: 40 to 180 min \pm 5% of display reading

Memory (with date and time): 100x2 users

Inflation: fuzzy logic controlled by electric pump

Deflation: automatic pressure release valve

Rating: 6VCC 4W

Power source: 4 "AA " batteries 1.5 V

Battery life (N° of measurements): 1000

Protection against electric shock: internally powered ME equipment (when using batteries)

Weight (console+cuff) without batteries: 460+163 g

Size: 191x120x85 mm

Cuff materials: nylon and polyester

