

OMRON INTELLI WRAP CUFF 22-42 cm HEM-FL31 for M3, M6, M7 Comfort - spare



Product Code	49920
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1233966
CND	Z12030285
EAN/UPC	4015672109056
GMDN	34978

Description

OMRON code: HEM-FL31-E
Accessories: Intelli adult wrap cuff
Size: 22-42 cm
Compatibility: M3 Comfort, M4 Intelli IT, M6 Comfort, M7 Intelli IT

Fits normal and large arms eliminating inaccurate readings associated with incorrect positioning.
AFib indication function.