



Product Code	49923
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	1774420
CND	Z1203020501
EAN/UPC	4015672111172

Description

OMRON RS7 INTELLI IT HEM-6232T-E - wrist

Intellisense technology makes measurement quick and comfortable. Advanced positioning sensor and cuff positioning guide for reliable measurements, which can be saved and tracked via Bluetooth through OMRON App.

Multilingual box and manual: GB, FR, IT, ES, DE, NL, RU, TR, Arabic.

Technical Specifications

Measuring method: oscillometric
 Measuring range:
 Pressure: 0 mmHg to 299 mmHg \pm 3 mmHg
 Pulse: 40 to 180 min \pm 5% of display reading
 Memory (with date and time): 100 \times 2 users
 Inflation: automatic
 Deflation: automatic
 Rating: DC 3 V 3.0 W
 Power source: 2 "AAA" batteries 1.5 V
 Battery life (No. of measurements): 300
 Protection against electric shock: internally powered ME equipment (when using batteries)
 Weight (console + cuff) without batteries: 91 g
 Size (mm): 91 \times 13 \times h 63
 Cuff materials: nylon and polyester
 Measurable circumference (cm): 13.5-21.5 (wrist)