

# ADAPTIC DIGIT 3M NON ADHERENT DRESSING - for finger - medium - sterile



Product Code	34994
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2301031
CND	M02030201
GTIN	25051978001437

## Description

3M™ ADAPTIC™ DRESSING - NON ADHERENT - sterile

- 3M Code: MAD013
- ADAPTIC™ DIGIT - for finger - medium
- Minimum order: box of 10

Non-adherent dressing that combines a low tack silicone coated material with a tubular bandage specifically designed for the treatment of finger and toe wounds.

Conformable and comfortable dressing, which allows for flexibility and maximum movement of the fingers and toes.

Easy to apply, helps minimize trauma associated with dressing changes.

It can be left on for several days.