

RI-CHAMPION SMART PRO+ BLOOD PRESSURE MONITOR WITH BLUETOOTH- 1735-BT - adult cuff



Product Code	49951
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2668630
CND	Z1203020501
EAN/UPC	4045396218895
PARAF	987475668

Description

RI-CHAMPION® SMART PRO+ BLOOD PRESSURE MONITOR WITH BLUETOOTH - 1735 BT - adult cuff

A clinically validated blood pressure monitor, designed for fast and reliable readings. Advanced performance with irregular heartbeat detection and averaging capabilities: analysis of 3 successive blood pressure measurements.

- Bluetooth® connectivity
- fast and quiet operation
- 200 measurements on one set of 4 AA batteries
- large, easy to read backlit screen
- meets EN1060-4, NIBP clinical investigation

Technical Specifications

Heart Rate Range: 40-199 beats/minute
Systolic Measurement Range: 60-255 mmHg
Diastolic Measurement Range: 30-195 mmHg
Pulse Rate Measurement Range: 40-199 beats / minute $\pm 4\%$ of reading
Maximum Inflation Pressure: 280 mmHg ± 3 mmHg or $\pm 2\%$ of reading
Power Source/Supply: 4x 1.5 V AA alkaline batteries / DC + 6 V / 1 A
Battery Life: 200 Readings
Memory: maximum 400 memory records
Size - weight: 141x121xh 72 mm - 350g without battery