



Product Code	50000
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

RDM (NSIS)	2683725
CND	Z1203020501
EAN/UPC	6970384252927

## Description

### MADRID ANEROID SPHYGMOMANOMETER

Professional, high quality sphygmomanometer with 0 to 300 mmHg scale.

Two handling modes for right and left-hand users.

Large dial Ø 59.5 mm with PVC shockproof ring.

Zinc-alloy handle with nickel plated.

Copper male and female connector.

Supplied with adult nylon cuff (22-32 cm) with D-ring, 1-tube PVC bladder, ribbed PVC bulb, black nylon bag and multilingual box and manual GB, FR, IT, ES, PT, DE, NL, SE, PL, CZ, HU, HR, RO, BG.

Bright dial

Male and Female connector

Left and Right hand use