

WUNDER RA300 PROFESSIONAL DIGITAL SCALE 300 kg - Class III - standard



Product Code	25003
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	2497351
CND	V0399
EAN/UPC	8023279250039
GMDN	38211

Description

WUNDER RA300 DIGITAL SCALE - 300 kg - Class III - standard

Professional digital column scale with a large and stable painted aluminium base with non-slip painted metal platform. The scale is equipped with a dual LCD display with triple simultaneous readout: weight, height and BMI. Adjustable display for clearer and more comfortable reading.
Functions: weight, tare, mother/child weight, hold, BMI, autozero, programmable auto-off.
4 large adjustable feet and 2 integrated wheels for easy transport.

Supplied with power adaptor, RS232 interface, level, 4 adjustable feet, 2 wheels and multilingual manual: GB, IT, ES. Available on gimaitaly.com: FR, PT, DE, NL, SE, NO, FI.
MDR certified.

Made in Italy.

Technical Specifications

WUNDER RA300 DIGITAL SCALE - 300 kg - Class III - standard
Capacity: 300 kg
Sensitivity: 100 g
Display: double LCD: LCD 20mm weight with 5 active digits; LCD 12 mm height and BMI
Scale size - weight: 355x570xh 950 mm - 9 kg
Packing size: 330x950x150 mm
Platform size: 355x310xh 75 mm
Column size: Ø60 x h 800mm
Power supply (included): adapter 12V - 2 A
Power supply (optional): rechargeable battery or 6xAAA alkaline batteries