

# WUNDER R150A PROFESSIONAL MECHANICAL COLUMN SCALE 150 kg - Class III



Product Code	25011
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	2516414
CND	V0399
EAN/UPC	8023279250114
GMDN	38267

## Description

### WUNDER R150A MECHANICAL COLUMN SCALE 150 kg - Class III

Professional mechanical column scale with sturdy and stable structure in strong aluminium coated with epoxy powder. Platform made of non-toxic, anti-slip rubber with removable mat for better sanitation.

Zero adjustment function.

Weight system with lever mechanism in cast aluminium and zinc plated steel knife edges to guarantee high performance over time. Large analogue dial Ø 180 mm, easy to read.

MDR certified.

Supplied with non-slip feet and user manual in GB, IT, ES. On request FR, DE.

Made in Italy.

## Technical Specifications

Capacity - sensitivity: 150 kg - 1000 g

Platform size: 240-260x300xh 85 mm

Column size: Ø 60xh 750 mm

Size - weight: 270x400xh 900 mm - 7 kg