

# OXY-0 PEDIATRIC FINGER OXIMETER - tiger - GB,FR,IT,PT,SE,RO



Product Code	35056
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2383371
CND	Z1203020408
EAN/UPC	8023279350562
PARAF	986897724

## Description

OXY-0 PEDIATRIC OXYMETER - tiger - GB, FR, IT, PT, SE, RO

Intended for measuring oxygen saturation (SpO<sub>2</sub>) and pulse rate of children at home, hospital, sports healthcare...  
LCD coloured display shows SpO<sub>2</sub> value, Pulse Rate value, Pulse Rate waveform and bar graph. 4 direction display with adjustable display mode and brightness. Automatic power off function within 5 seconds.  
It includes one hanging rope, two AAA batteries and user manual.

## Technical Specifications

Display: LCD coloured display

SpO<sub>2</sub>:

- range: 0% ~ 100%, (resolution is 1%)
- accuracy: 70% ~ 100%: ±2%, < 70% unspecified

Pulse Rate

- range: 30 bpm ~ 250 bpm (resolution 1 bpm)
- accuracy: ±2 bpm or ±2% select larger

Power: 2xAAA 1.5 V alkaline batteries

Battery working hour: two batteries can work continuously for 20 hours

Dimension - weight: 59x36xh 35 mm - 50g (including batteries)