



Product Code	35066
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2060328
CND	Z1203020408
EAN/UPC	8023279350661
UMDNS	17148
PARAF	986897700
GMDN	45607

Description

O2ring™ continuous monitoring oximeter - kid

Innovative, wearable, adaptable and reliable pulse oximeter, designed to monitor heart rate, oxygen saturation and perfusion index.

Monitor oxygen level anytime, anywhere

- comfortable wearing experience
- track and record oxygen level and heart rate
- smart vibration triggered by low SpO2
- built-in rechargeable battery for overnight use
- sync real-time data with App via Bluetooth

Track and record oxygen level overnight

During your sleep, once a low oxygen level or an abnormal heart rate is detected, O2Ring will vibrate immediately and softly to remind you. After syncing the recording data to App, detailed analysis and trends of your oxygen level and heart rate can be reviewed.

Supplied with data/charging cable, user manual (GB, FR, IT , ES , DE).

PC software can be downloaded online.

Technical Specifications

Blood oxygen: range: 70-99% / accuracy: 80-99%: $\pm 2\%$, 70-79%: $\pm 3\%$

Pulse rate: range: 30 to 250 bpm - ± 2 bpm or $\pm 2\%$

Vibration source: Low oxygen level, high/low pulse rate

Recorded parameters: Blood oxygen level, pulse rate, exercise

Weight: 15 g

Size: device 38x30x38 mm; ring 50-82 mm

Built-in memory: 4 sections, up to 10 hours each section

Battery: Rechargeable lithium battery 3.7 Vdc

Use time: 12-16 hours

Charge time: 2-3 hours

Bluetooth: Bluetooth 4.0 BLE

Compatibility: iOS 9.0+; Android 5.0+

Protection: electric shock: Type BF electromagnetic: Group I, Class B dust and waterproof rating: IP24