

# NONIN ONYX VANTAGE PULSE OXIMETER - 9590 - black



|               |                |
|---------------|----------------|
| Product Code  | 35085          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II B           |

|            |               |
|------------|---------------|
| RDM (NSIS) | 451807        |
| CND        | Z1203020408   |
| EAN/UPC    | 0833166004878 |
| PARAF      | 987475504     |
| GMDN       | 17148         |

## Description

### NONIN ONYX® VANTAGE 9590 PULSE OXIMETER - black

The only finger pulse oximeters with scientifically proven accuracy in motion and low perfusion. Tested for use on fingers, thumbs, or toes.

- Accurate - Scientifically proven performance in the widest range of patient populations and settings.
- Durable - Protects against dropping and water spills; allows for thousands of uses in any environments and exceed IP32 water ingress testing.
- Efficient - Up to 6000 spot checks or 36 hours of continual operation on two AAA batteries.
- Safe - Built in the USA with lead-free and latex-free construction.
- Versatile - suitable for pediatric to larger adult patients. Wide range of digit thicknesses from 8 mm to 25.4 mm.
- Warranted - Industry-leading 4-year warranty.
- Cost Effective - Proven accuracy, quality and durability make the Onyx Vantage a smart purchase for the long term.

Multilingual CD manual: GB, FR, IT, DE, ES, PT, GR, NL, PL, NO, FI, DK, SE.

Large bright LED Display. SpO2 and pulse rate can be read from any angle, day or night.  
Lead (PB) free and latex free

## Technical Specifications

Pulse rate range: 18-321 bpm  $\pm 3$  digit (20-250 bpm)  
SpO2 range: 0-100%  $\pm 2$  digit (70-100%)  
Low perfusion SpO2: 70-100%  $\pm 2$  digit  
Low perfusion pulse rate: 40-240 bpm  $\pm 3$  digit