

NONIN PALM SAT 2500A PULSE OXIMETER



Product Code	35086
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	287228
CND	Z1203020408
EAN/UPC	0849686090178
PARAF	962357099
GMDN	17148

Description

PALM SAT® 2500A PULSE OXIMETER with rubber reusable adult probe and alarm

A compact versatile and reliable pulse oximeter designed and tested to exceed all durability requirements.

- Pure Sat® technology offers accuracy & speed in getting a reading, even with varying skin pigmentations or polish and is able to monitor, with declared accuracy in challenging conditions, such as motion & low perfusion, for all patients
- it is tested to assure performance in challenging conditions, such as drop and liquid ingress (IPX 2). Granted by US Army/Air Force Aeromedical certification

Audible and visual alarms, 72 hours memory, long battery life (60 hours) or optional NiMH rechargeable battery for 40 hours.

- Bright LED display of SpO2 and pulse rate
- Three-coloured signal quality indicator
- Blinking alarm silence indicator Low battery indicator
- Rugged rubberized case with IPX 2 level of liquid ingress protection

Size: 7x13.8x3.2 cm
Weight: 213 g

Supplied with multilingual cd-rom manual: GB, FR, IT, DE, ES, PT, GR, NL, PL, NO, FI, DK, SE.

Made in the USA.

Technical Specifications

SpO2 range: 0% to 100% ±2 (70-100% without movement)
Pulse rate range: 18 to 321 bpm ±3 (18-300% without movement)
Power: 4AA alkaline batteries 1.5 V (for 60 hours use)