

# COLOURED TRAD STETHOSCOPE - blue



Product Code	51002
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2513131
CND	C9005
EAN/UPC	8023279510027
PARAF	977817511
GTIN	08023279510027
GMDN	13755

## Description

COLOURED TRAD STETHOSCOPES - LATEX FREE

- TRAD single head - blue

Coloured aluminium stethoscopes, lightweight and reliable.  
Chestpiece diaphragm Ø 43.5 mm, bell Ø 30 mm external sizes.

Brass binaural.

Tubing: length 62 cm, int. Ø 4 mm, ext. Ø 8 mm.

Supplied in cardboard box.

Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic